

CYCLIC AND LOAD TO FAILURE PROPERTIES OF ALL-SUTURE ANCHORS IN SYNTHETIC ACETABULAR CANCELLOUS BONE

Nathan P. Douglass, MD
Anthony W. Behn, MS
Marc R. Safran, MD

Department of Orthopedic Surgery



STANFORD
UNIVERSITY

ISHA Annual Meeting
San Francisco, USA
September 2016


Disclosures

- Dr. Safran
 - Consultant - ConMed and Medacta
 - Fellowship support – Smith & Nephew, Ossur, and ConMed Linvatec
 - Research Support - ConMed Linvatec, ISAKOS, Ferring Pharmaceuticals
 - Royalties – Smith & Nephew, Stryker, DJO
- Dr. Douglass and Mr. Behn have no disclosures




Introduction

- Metal, polyether ether ketone (PEEK) and various biocomposites have been routinely used in shoulder labral repair
- However, there is a paucity of data on the performance of the more recently-developed all-suture anchors (ASAs) [1-3]
- In particular, ASAs have not been tested in acetabular bone, or bone substitute simulating the range of bone properties found in the acetabulum.



Questions

- 1. Do all-suture anchors have better or worse fixation than a commonly used standard suture anchor?
- 2. Is there variation in fixation between all-suture anchors of different designs?
- 3. Does high density bone affect the fixation of these various suture anchors?



Methods

Test all ASAs marketed for hip labral repair

- 12 anchors
- Bioraptor 2.3 PK = control (PEEK ribbed implant)
- 7-11 anchors per group based on previous studies [1-3]

Sawbones with properties replicating:

- High density – acetabular rim (30 pcf)
- Low density – glenoid (20 pcf)
- Based on multiple cadaveric studies [4-12]


Comparing outcomes in 20 pcf vs. 30 pcf for each anchor

- Welch T-test

Comparing anchors to each other, in both densities

- Welch ANOVA with Games-Howell post hoc test


Comparisons only made when at least 3 anchors survived (or were present)



Commercially-Available All-Suture Anchors


Name	Company	Drill size (mm)	Suture size	# suture strands	Indications
Suturefix Ultra 1.7	Smith & Nephew	1.7	2	1	hip and shoulder labrum
Suturefix Ultra 1.9	Smith & Nephew	1.9	1	2	shoulder labrum
JuggerKnot 1.0	Biomet Sports Medicine	1.0	2-0 or 3-0	1	hand, small joint
JuggerKnot 1.4	Biomet Sports Medicine	1.4	1	1	shoulder labrum
JuggerKnot 1.45 #1	Biomet Sports Medicine	1.45	1	1	hip labrum
JuggerKnot 1.45 #2	Biomet Sports Medicine	1.45	2	1	hip labrum
JuggerKnot 1.5	Biomet Sports Medicine	1.5	2	1	shoulder labrum
JuggerKnot 2.1	Biomet Sports Medicine	2.1	1	1	shoulder labrum
JuggerKnot 2.9	Biomet Sports Medicine	2.9	2	2	rotator cuff
Y-Knot 1.3	ConMed Linvatec	1.3	2	1	hip and shoulder labrum
Y-Knot 1.8	ConMed Linvatec	1.8	2	2	shoulder labrum
Iconix 1	Stryker Endoscopy	1.4	2	1	broad, including hip
Iconix 2	Stryker Endoscopy	2.3	2	2	broad, including hip
Iconix 3	Stryker Endoscopy	2.3	2	3	broad, including hip
Iconix 25	Stryker Endoscopy	2.3	5	2	broad, including hip
Q-Fix 1.8	ArthroCare	1.8	2	1	hip and shoulder labrum
Q-Fix 2.8	ArthroCare	2.8	2	2	rotator cuff

Red arrowhead indicates anchor tested in this study



Methods


- Deployment Force: 40 N
- Cyclic loading protocol stages
 1. 10 – 50 N for 200 cycles, 0.5 Hz
 2. 10 – 100 N for 200 cycles
 3. Load to failure, 10 mm/s



Instron testing setup




Peak Displacement Cycle 200



- In 20 pcf test blocks, the Bioraptor 2.3 PK and Q-Fix 1.8 showed significantly ($p < 0.002$) less peak displacement than all other anchors
- In 30 pcf test blocks, the Q-Fix 1.8 showed a trend towards less displacement than the Bioraptor 2.3 PK ($p = 0.11$) but not Iconix 25, and significantly less than all others ($p < 0.02$).




Post-Cyclic Displacement




- In 20 pcf test blocks, the Q-Fix 1.8 displaced less ($p < 0.01$) than all other anchors.
- In 30 pcf test blocks, the Q-Fix 1.8 displaced less than all anchors ($p < 0.02$), except the Iconix 25 ($p = 0.18$)



Maximum Load




- Higher failure loads in 30 pcf than 20 pcf for all anchors, except Iconix 2 ($p = 0.09$)
- No statistically significant differences among top performers in either 20 or 30 pcf blocks
- Anchors routinely failed by suture pullout, except the Q-Fix 1.8 which failed by suture breakage in 30 pcf blocks




Limitations

- Synthetic bone substitute
- Non-aqueous environment
- Absence of cortical layer
 - Does not affect majority of rotator cuff anchors [1, 13, 14]
 - However, cortical layer often disrupted by acetabuloplasty or preparation
- 400 Cycles
 - Based on previous studies Barber and Herbert
 - Hip / Shoulder often braced post-operatively
- Unknown clinically relevant load to failure, displacement
 - 250N in rotator cuff [15-18]



Discussion

- ASAs have improved *fixation* in higher density bone
 - Did not test *deployment* directly with imaging
- ASAs generally performed worse than control standard anchor
 - With exceptions - not in 20 pcf, with over 50 N load
- Q-Fix 1.8 outperformed all anchors including the control anchor in all testing conditions



References

1. Barber FA, Herbert MA. Cyclic loading biomechanical analysis of the pullout strengths of rotator cuff and glenoid anchors: 2013 update. *Arthroscopy*. 2013;29(5):832-844.
2. Mazzocca AD, Chowanec D, Cote MP, Fierra J, Apostolakos J, Nowak M, et al. Biomechanical evaluation of classic solid and novel all-soft suture anchors for glenoid labral repair. *Arthroscopy*. 2012;28(5):642-648.
3. Galland A, Airaudi S, Gravier R, Le Cann S, Chabrand P, Argenson J-N. Pullout strength of all suture anchors in the repair of rotator cuff tears: a biomechanical study. *Int Orthop*. 2013;37(10):2017-2023.
4. Dalstra M, Huiskes R, Odgaard A, van Erning L. Mechanical and textural properties of pelvic trabecular bone. *J Biomech*. 1993;26(4-5):523-535.
5. Helgason B, Perilli E, Schileo E, Taddei F, Brynjólfsson S, Viceconti M. Mathematical relationships between bone density and mechanical properties: a literature review. *Clin Biomech (Bristol, Avon)*. 2008;23(2):135-146.
6. Speirs AD, Beaulé PE, Rakhra KS, Schweitzer ME, Frei H. Increased acetabular subchondral bone density is associated with cam-type femoroacetabular impingement. *Osteoarthr Cartil*. 2013;21(4):551-558.
7. Kalouche I, Crépin J, Abdelmoumen S, Mitton D, Guillot G, Gagey O. Mechanical properties of glenoid cancellous bone. *Clin Biomech*. 2010;25(4):292-298. doi:10.1016/j.clinbiomech.2009.12.009.
8. Mansat P, Barea C, Hobatho MC, Darmana R, Mansat M. Anatomic variation of the mechanical properties of the glenoid. *J Shoulder Elbow Surg*. 1998;7(2):109-115.
9. Frich LH, Jensen NC, Odgaard A, Pedersen CM, Sjøbjerg JO, Dalstra M. Bone strength and material properties of the glenoid. *J Shoulder Elbow Surg*. 1997;6(2):97-104.
10. Frich LH, Odgaard A, Dalstra M. Glenoid bone architecture. *J Shoulder Elbow Surg*. 1998;7(4):356-361.
11. Adler E, Stuchin SA, Kummer FJ. Stability of press-fit acetabular cups. *J Arthroplasty*. 1992;7(3):295-301.
12. Bardana DD, Burks RT, West JR, Greis PE. The effect of suture anchor design and orientation on suture: An in vitro study. *Arthroscopy*. 2003;19(3):274-281.
13. Barber FA, Herbert M a., Beavis RC, Barrera Oro F. Suture anchor materials, eyelets, and designs: update 2008. *Arthroscopy*. 2008;24(8):859-867.
14. Barber AF, Herbert MA, Hapa O, Rapley JH, Barber CAK, Bynum JA, et al. Biomechanical analysis of pullout strengths of rotator cuff and glenoid anchors: 2011 update. *Arthrosc - J Arthrosc Relat Surg*. 2011;27(7):895-905.
15. Burkhart SS, Diaz Pagán JL, Wirth MA, Athanasiou KA. Cyclic loading of anchor-based rotator cuff repairs: confirmation of the tension overload phenomenon and comparison of suture anchor fixation with transosseous fixation. *Arthroscopy*. 1997;13(6):720-724.
16. Schneeberger AG, von Roll A, Kalberer F, Jacob HAC, Gerber C. Mechanical strength of arthroscopic rotator cuff repair techniques: an in vitro study. *J Bone Joint Surg Am*. 2002;84-A(12):2152-2160.
17. Gerber C, Schneeberger AG, Beck M, Schlegel U. Mechanical strength of repairs of the rotator cuff. *J Bone Joint Surg Br*. 1994;76(3):371-380.
18. Mazzocca AD, Millett PJ, Guancho CA, Santangelo SA, Arciero RA. Arthroscopic single-row versus double-row suture anchor rotator cuff repair. *Am J Sports Med*. 2005;33(12):1861-1868.

